

Western Blot

RR-WB#

Samples:

Purpose:

Making the gels:

- Make up separating gel:

Volume for Different Percentages of Acrylamide					
Component	6%	8%	10%	12%	15%
lower gel 4X buffer	2.5ml	2.5ml	2.5ml	2.5ml	2.5ml
dH ₂ O	5.9ml	5.4ml	4.9ml	4.4ml	3.65ml
40% Acrylamide	1.5ml	2.0ml	2.5ml	3.0ml	3.75ml
10% APS*	100.0µl	100.0µl	100.0µl	100.0µl	100.0µl
TEMED*	10.0µl	10.0µl	10.0µl	10.0µl	10.0µl
kDa Range:	60-200		16-70		12-45

*Add just prior to pouring the gel

- Pour into gel apparatus using a P1000 and quickly overlay with 1 mL of hydrated Butanol
(Leave $\frac{1}{4}$ space for stacking gel.)
- Let sit for 30-45 min
- After polymerization, remove the Butanol overlay via vacuum or by gravity
- Rinse with dH₂O
- Make up stacking gel:

4% Stacking gel	
Component	Volume
upper gel 4X buffer	2.5ml
dH ₂ O	6.4ml
40% Acrylamide	1.0ml
10% APS*	100µl
TEMED*	10µl

*Add just prior to pouring the gel

- Insert the comb and let sit to polymerize for at least 30 min
- Place the gel in the electrophoresis chamber, cover with running buffer, remove the comb and rinse the wells with buffer to remove unpolymerized acrylamide


Western Blot

Preparing the samples:

- Add Laemmli buffer to aliquoted samples
- Incubate the sample at 95°C for 5 minutes to denature the proteins
- Mix by vortexing, blip spin and load onto the gel
- Add sample buffer into all empty wells for an even dye front
- Run gel at 180V until the dye front has emerged from the gel (about 1 hour)

Blotting:

- Remove stacking layer from the gel and clip the gel on the top, left corner by the marker
- Place membrane in tray filled with dH₂O
- Fill large tray with transfer buffer and assemble the blot:



**place the gel down so that the cut is on the top, RIGHT corner*

**place the membrane so that the shiney side is up and the matte side is down on the gel*

**use fingers to "walk" over gel, but do not rub*

- Flip the blot assembly over and place in semi-dry blotter (membrane-side down)
- Set to 100-300 mA, 15 V and 2-5 Watts (with a CONSTANT amp):

2 mini gels = 80 minutes

4 mini gels = 120 minutes

6 mini gels = 240 minutes

- Stain with 0.1% Ponceau S (in 1% acetic acid) in order to visualize protein bands by swirling the solution onto the blot for 5 minutes
- Rinse blot with dH₂O to remove background staining
- Scan or photograph the blot to confirm equal loading in each lane

**Ponceau S will come off during the blocking step.*

***Blot may be rinsed in dH₂O and stored at RT *prior* to Ponceau stain.*

****Blot may be blocked, washed with 1X PBS and then stored at 4°C. Use dH₂O to wash if using an AP-conjugated secondary antibody.*

Blocking:

- Block in 5% non-fat dry milk in TBST at RT for 30 min or at 4°C O.N.

Western Blot

Label Blots:

